


SILVER CAP (-U.S.A.-) Glass-body pH electrodes:



Common features of SILVER CAP (-U.S.A.-) Glass-body pH electrodes:

- Ag/AgCl reference electrode if applicable
- Annular Ceramic Junction except electrode F which is wick.
- Temperature Range: 0 - 65 °C(Refillable) & 0 – 45 °C (Sealed)
- BNC connector with 3 ft (~ 1 meter) low-noise cable. Other connectors are available, please specify.

Note: The bulb color and shape may differ from the pictures.

P/N	Figure	Half / Combo	pH range	Bulb Shape	Sealed / Refillable	Junction type	Dimensions mm
PH5001-3B	E	Combo	0-14	Ball	Sealed	Single	12 x 120
PH5005-3B	E	Combo	0-14	Ball	Refillable	Single	12 x 120
PH5007-3B	E	Combo	0-14	Dome	Refillable	Single	12 x 120
PH5009-3B	E	Combo	0-14	Dome	Sealed	Single	12 x 120
PH5012-3B	H	Combo	0-14	Ball	Sealed	Double	12 x 120
PH5016-3B	H	Combo	0-14	Ball	Refillable	Double	12 x 120
PH5020-3B	A	Combo	0-14	Ball	Refillable	Double	8.2x60 semi portion
PH5025-3B	A	Combo	0-14	Ball	Sealed	Single	8.2x60 semi portion
PH5027-3B	A	Combo	0-14	Ball	Refillable	Single	8.2x60 semi portion
PH5028-3B	A	Combo	0-14	Ball	Sealed	Double	8.2x60 semi portion
PH5030-3B	B	Combo	0-14	Cylinder	Refillable	Double	8.2x65 semi portion
PH5035-3B	B	Combo	0-14	Cylinder	Sealed	Single	8.2x65 semi portion
PH5037-3B	B	Combo	0-14	Cylinder	Refillable	Single	8.2x65 semi portion
PH5038-3B	B	Combo	0-14	Cylinder	Sealed	Double	8.2x65 semi portion
PH5040-3B	C	Combo	0-14	Spear	Refillable	Double	8.2x65 semi portion
PH5045-3B	C	Combo	0-14	Spear	Sealed	Single	8.2x65 semi portion

PH5047-3B	C	Combo	0-14	Spear	Refillable	Single	8.2x65 semi portion
PH5048-3B	C	Combo	0-14	Spear	Sealed	Double	8.2x65 semi portion
PH5050-3B	D	Combo	0-14	Ball	Refillable	Double	8.2 x110 semi portion
PH5057-3B	D	Combo	0-14	Ball	Refillable	Single	8.2 x110 semi portion
PH5065-3B	F	Combo	0-14	Dome	Refillable	Single	12 x 120
PH5067-3B	F	Combo	0-14	Ball	Refillable	Single	12 x 120
PH5070-3B	G	Combo	0-14	Ball	Refillable	Double	12 x 120
PH5077-3B	G	Combo	0-14	Ball	Refillable	Single	12 x 120
PH5081-3B	I	Half	0-14	Ball	N/A	N/A	12 x 120
PH5083-3B	J	Half	0-14	Dome	N/A	N/A	12 x 120