

ALTRONIX TPA-V

pH/mV/Temp Meter



Model: TPA-V

Features:

1. Microprocessor based with splash proof housing.
2. Simultaneous display pH and Temperature.
3. Automatic or manual temperature compensation. °C/°F switchable
4. Simple to calibrate by one keyboard for 3 points buffer.
5. Calibration value can be adjusted as needed.
6. Indicate percentage of slope (PTS) after calibration.
7. Low battery and consumption indicator.
8. Auto shut off after 10 minutes of non use.
9. Rubber protective holster with magnetic.

Specifications:

	TPA-V			
		pH	mV	Temp.
Range:		0-14.00	±1999	0 -100°C
Accuracy:		±0.01+1 digit	±1+1 digit	±0.2°C+1digit
Resolution:		0.01	1	0.1°C
Compensation:		ATC 0- 100°C		
Calibration:	pH 4.00, 7.00, 10.00			
Battery:	9V			
Dimensions:	108 x 75 x 30 mm			
Weight:	135 g (with battery)			

Accessories for complete kit:

1. 9V Battery.
2. Instruction manual.
3. Gift box.
4. pH electrode with Pt 100 (Optional).
5. pH electrode (U.S.A.).
6. ATC(Optional)
7. ORP electrode (U.S.A.) (Optional).
- 8.- Rubber protective holster with magnetic.