

Digital display Orbital Shaker Technical Data

Picture						
Model	OS-20F	OS-40F	LS-40F	OS-20Pro	OS-40Max	LS-40Max
Motor type	Brushless motor					
Plate Size	250*250mm	320*320mm	320*320mm	250*250mm	320*320mm	320*320mm
Turn	Circular	Circular	Linear and reciprocating	Circular	Circular	Linear and reciprocating
Max Load	2kg	3kg	3kg	2kg	5kg	5kg
Orbit Diameter	20mm	20mm	Linear and reciprocating 20mm	20mm	20mm	Linear and reciprocating 20mm
Speed Range	40-200rpm					
Timing range	1min - 23h59min					
Voltage	100-240V, 50/60Hz					
Allowable Ambient Temp	5~40°C 80%RH					
Size	355*260*100mm	435*330*140mm	435*330*140mm	355*260*135mm	435*330*140mm	435*330*140mm
Net weight	3.65kg	3.75kg	4kg	4.85kg	8.1kg	8.1kg