

## REFRACTOMETROS LANCET

**Lancet**  
INSTRUMENTS

### OPTICOS

Itemº	ESCALA	RANGO	DIVISIÓN	PRECISIÓN
<b>RHB-10 ATC</b>	BRIX	0/10%	0,10%	+/-10%
<b>RHB-18 ATC</b>	BRIX	0/18%	0,10%	+/-10%
<b>RHB-32 ATC</b>	BRIX	0/32%	0,20%	+/-10%
<b>RHB-62 ATC</b>	BRIX	28/62%	0,20%	+/-10%
<b>RHB-82 ATC</b>	BRIX	45/82%	0,50%	+/-50%
<b>RHB-90 ATC</b>	BRIX	58/90%	0,50%	+/-50%
miel	BAUME	38-43 Be	0,5 Be	+/-0,50Be
	AGUA	12/27%	1%	+/-1%
<b>Ref 107</b>	BRIX	0/90% 3 Escalas	0,20%	+/-0,2%
<b>Ref 108</b>	BRIX	0/80% 2 Escalas	1%	+/-1%
<b>Ref 108</b>	BRIX	0/80% 1 Escala	1%	+/-1%
<b>RHS-10 ATC</b>	SALINA	0-10%	0,10%	+/-0,10%1(ppt)
<b>RHS-28 ATC</b>	SALINA	0-28%	0,20%	+/-0,2%2(ppt)
<b>RHC-200 ATC</b>	SP:	0/12 g/dl	0,2 g/dl	+/-0,2 g/dl ATC
Clinico	SG:	1000/1050sg	0,005 sg	+/-0,005sg
	RI:	1,3330/1,3600 Ri	0,00005RI	+/-0,0005RI
<b>RHW-25 ATC</b>	Alcohol	0-25% VOL	0,20% VOL	0,2% VOL
<b>RHW-25 ATC</b>	Alcohol	0/25%	0,20% VOL	0,2% BRIX
	BRIX	0/40% BRIX	0/40% BRIX	0,2% VOL
<b>RHA-200 ATC</b>	BATERIA	E:-60°C-0°C	5°C%	5°C +/-
bateria/fluidos		P:-50°C-0°C	5°C%	5°C +/-
		B:1.1000-1400sg	0,01sg%	0,01sg +/-
<b>Ref 601</b>	Gemas			
<b>Opcional</b>	Con luz de led más:			
<b>2 WJA</b>	Rango 1300/1700 nD	Resolucion 0,0005 Nd	ABBE de Mesa	
	Rango 0/95%	Resolucion 0,25%		

### DIGITALES

Itemº	ESCALA	RANGO	DIVISIÓN	PRECISIÓN
<b>LARB-45</b>	BRIX	0-45%	0,1%	+/-0,1%
	IND.REF	1,3330-1,4098	0,0001	+/-0,0003
<b>LARB-28-65</b>	BRIX	28-65%	0,1%	+/-0,1%
	IND.REF	1,3800-1,4535	0,0001	+/-0,0003
<b>LARB-60-92</b>	BRIX	60/92%	0,1%	+/-0,1%
	IND.REF	1,4400-1,5230	0,0001	+/-0,0003
<b>LARS-95</b>	BRIX	0-95%	0,1%	+/-0,1%
	IND.REF	1,4400-1,5230	0,0001	+/-0,0003
<b>LARS-28</b>	SALINA	0-28%	0,1%	+/-0,1%
<b>SALINO</b>	IND.REF	1,3300-1,3900	0,0001	+/-0,0003
<b>LARC-300</b>	Orina SP.G	1,000-1,050	0,001	+/-0,001
<b>CLINICO</b>	Serum P.	0-12g/dL	0,1	+/-0,1
	IND.REF	1,3330-1,3900	0,0001	+/-0,0003
<b>LARW-400</b>	BRIX	0,35%	0,1%	+/-0,1%
<b>VINO</b>	IND.REF	1,4400-1,5230	0,0001	+/-0,0003