

pH electrodes



PY41



GL42



PS45



TP-46T30K



PF47



ID4410



ID4520



ID4520F



Model	PY-41	GL-42	PS-45	TP-46T30K	PF-47	ID-4410	ID-4520	ID-4520F
Description	Plastic electrode	Glass electrode	Spear electrode	With 30K temp. sensor	Flat electrode	on-line	on-line	on-line/flat
Range	0-14	0-14	0-14	0-14	0-14	0-14	0-14	0-14
Temp.	0-80°C	0-90°C	0-80°C	0-80°C	0-80°C	0-90°C	0-80°C	0-80°C
Type Ref	Ag/AgCl	Ag/AgCl	Ag/AgCl	Ag/AgCl	Ag/AgCl	Ag/AgCl	Ag/AgCl	Ag/AgCl
Junction	Fiber	Ceramic	PVC	Ceramic	Cricoid Teflon	Ground glass sleeve	Cricoid Teflon	Cricoid Teflon
Refillable	No	Yes	No	No	No	Yes	No	No
Body	PC	Glass	PC	PC	PC	Glass	GFRPP	GFRPP
Connectr	BNC	BNC	BNC	BNC + 3.5 ψPlug or RCA	BNC	Spade Lug or BNC	Spade Lug or BNC	Spade Lug or BNC
Size(mm)	160 x 12	160 x 12	174 x 12	160 x 12	160 x 12	160 x 12	176 x 16	176 x 16
Cable	1 m	1 m	1 m	1 m	1 m	1.2 m	5 m	5 m