

REFRACTOMETROS LANCET

Lancet
INSTRUMENTS

OPTICOS

Itemº	ESCALA	RANGO	DIVISIÓN	PRECISIÓN
RHB-10 ATC	BRIX	0/10%	0,10%	+/-10%
RHB-18 ATC	BRIX	0/18%	0,10%	+/-10%
RHB-32 ATC	BRIX	0/32%	0,20%	+/-10%
RHB-62 ATC	BRIX	28/62%	0,20%	+/-10%
RHB-82 ATC	BRIX	45/82%	0,50%	+/-50%
RHB-90 ATC	BRIX	58/90%	0,50%	+/-50%
miel	BAUME	38-43 Be	0,5 Be	+/-0,50Be
	AGUA	12/27%	1%	+/-1%
Ref 107	BRIX	0/90% 3 Escalas	0,20%	+/-0,2%
Ref 108	BRIX	0/80% 2 Escalas	1%	+/-1%
Ref 108	BRIX	0/80% 1 Escala	1%	+/-1%
RHS-10 ATC	SALINA	0-10%	0,10%	+/-0,10%1 (ppt)
RHS-28 ATC	SALINA	0-28%	0,20%	+/-0,2%2 (ppt)
RHC-200 ATC	SP:	0/12 g/dl	0,2 g/dl	+/-0,2 g/dl ATC
Clinico	SG:	1000/1050sg	0,005 sg	+/-0,005sg
	RI:	1,3330/1,3600 Ri	0,00005RI	+/-0,0005RI
RHW-25 ATC	Alcohol	0-25% VOL	0,20% VOL	0,2% VOL
RHW-25 ATC	Alcohol	0/25%	0,20% VOL	0,2% BRIX
	BRIX	0/40% BRIX	0/40% BRIX	0,2% VOL
RHA-200 ATC	BATERIA	E:-60°C-0°C	5°C%	5°C +/-
bateria/fluidos		P:-50°C-0°C	5°C%	5°C +/-
		B:1.1000-1400sg	0,01sg%	0,01sg +/-
Ref 601	Gemas			
Opcional		Con luz de led más:		
2 WJA	Rango 1300/1700 nD	Resolucion 0,0005 Nd	ABBE de Mesa	
	Rango 0/95%	Resolucion 0,25%		

DIGITALES

Itemº	ESCALA	RANGO	DIVISIÓN	PRECISIÓN
LARB-45	BRIX	0-45%	0,1%	+/-0,1%
	IND.REF	1,3330-1,4098	0,0001	+/-0,0003
LARB-28-65	BRIX	28-65%	0,1%	+/-0,1%
	IND.REF	1,3800-1,4535	0,0001	+/-0,0003
LARB-60-92	BRIX	60/92%	0,1%	+/-0,1%
	IND.REF	1,4400-1,5230	0,0001	+/-0,0003
LARS-95	BRIX	0-95%	0,1%	+/-0,1%
	IND.REF	1,4400-1,5230	0,0001	+/-0,0003
LARS-28	SALINA	0-28%	0,1%	+/-0,1%
SALINO	IND.REF	1,3300-1,3900	0,0001	+/-0,0003
LARC-300	Orina SP.G	1,000-1,050	0,001	+/-0,001
CLINICO	Serum P.	0-12g/dL	0,1	+/-0,1
	IND.REF	1,3330-1,3900	0,0001	+/-0,0003
LARW-400	BRIX	0,35%	0,1%	+/-0,1%
VINO	IND.REF	1,4400-1,5230	0,0001	+/-0,0003